

A complete prepared system for multi organ preparation

General information about products from Meddev

- A specific designed system
- All parts can be connected fast and easy
- All products are made of PVC. No latex, which can lead to allergic reactions for some patients.
- High quality products. Produced in Switzerland according to European standards.
- Most products are double packed, for extra safety in the operating room.
- Custom made products for your hospital is available on request
- Fast delivery out of Nordic stock placed in Copenhagen.



Organ Perfusion Sets: Double line and Single line

For infusion of preservation solutions. The length enables for space among transplant teams during harvesting of organs. Allows for fast flow of preservation solution. Spikes allows for quick and easy connection to preservation solution bag. The drip chamber is large for easy viewing and control. Air filter prevents any air bubbles to pass on to donate organ. Roller clamps allows for controlled flow of solution fast, medium and slow. End connector allows for connection to any type of cannula.

Custom sets or individual modifications to the standard perfusion set are available. Minimum order requirement for special manufactured products is 100 units.



Explant cannulas:

The cannula leads the perfusion solution into the artery. It is connected directly to the perfusion set. It is also used for draining the organ for fluids. Available in pediatric and adult diameters Fr.8 - Fr.32. Insertion is easy into vessel or artery. The hub allows for the surgeon to tie suture tightly, allowing the cannula to remain in position at the point of cannulation. Other cannulas may “pop” out due to pressure. The point will not damage the vessel or artery. Cannula end can be connected to any tubing.



Back Table Flush adapter

A small cannula for quick and easy cannulation of the renal artery. Easily attached to any type of tubing which is then connected to the perfusion set.



Portal cannula

Portal cannula for quick and easy perfusion of the portal vein of the liver.
Available in 3 sizes. Fr. 7, Fr. 10 and Fr. 12.

Vena cava Drainset:

A disposal set that helps to maintain a clean operating room floor during organ donation. All fluids are collected into the bag and then disposed neatly. Prevents overflow. 5 liter or 10 liter bags allows for a high volume of fluids. Flow is due to gravity, no pump needed. Explant cannula Fr. 32, Fr. 34 or Fr. 36 is available for connection to donor trough the vena cava.



Cardioplegia set with or without needles

Same principles apply due to organ perfusionsets. Sets are available with needles for ventilation of bottles. Sets without needles are available for plastic bags.